

GeneTel Laboratories, LLC *Leader in Proteomics Abs and Abs Arrays*
1202 Ann Street, Madison, WI 53713 Phone 877-248-4316 Fax 608-257-1385 Email info@genetel-lab.com

Product: Biotinylated Chicken Anti-H3 Relaxin B chain from human, affinity purified.

Catalog Number: OB0046

Lot Number: 130209

Size: 100ul at 0.5mg/ml in PBS with 0.02% Na Azide


Immunogen: Human Relaxin H3 B chain peptide KWGCSKSEISSLC.

Specification: OB0046 specifically binds to human H3 Relaxin B chain. The product may cross-react with H3 relaxin of other species.

Suggested Usage: The product can be used at 1:2,000 to 1:10,000 for ELISA tests in combination with OB0044 anti-H3 A chain. User needs to optimize the concentrations in assays. The product can be used in combination with GeneTel's streptavidin conjugates such as SAHRP (Streptavidin HRP) and SAAP (Streptavidin AP).

Short Term Storage: 2-8°C for one week.

Long Term Storage: -20°C as an undiluted liquid for up to two year.
Please feel free to call us if you have any question regarding to this product or any other GeneTel's product.



H3 ELISA Test using OB0044 as coating Abs and OB0046 as detection Abs